

MTP-1000 Automatic Tissue Processor



Features:

1. Chinese menu, man-machine interaction, simple operation.
2. Prompt judge and automatic processing in case of abnormal condition, which can prevent the tissue from being damaged.
3. Double-control heat without water
4. Based on different tissues, several running programs available for customers
5. Self-provided power supply can ensure to continue 6 hours' working even if the electricity has been cut off.
6. Environment friendly: The air as purified through internal circulation can be exhausted off the window.
7. With transmission of gear and rack, mechanical failures can be reduced.
8. Several standby programs available

Technical Data

9. 12 tanks (3 wax cylinders in the end)
10. Capacity of single- tanks: 1000ml
11. Basket can hold 60 specimens
12. Single tanks processing time: set arbitrarily in the range of 0 to 99 hour(s)

Accessories:

13. Standard Setting List

- 14. Reagent bottle 9 pcs
- 15. 2. Tissue basket 1 pc
- 16. 3.The exhaust pipe 1 pc
- 17. 4.Tissue cassette 20 pcs
- 18. 5. Fuse (5 A) 2 pcs
- 19. 6. Power line 1 line
- 20. Optional Accessories
- 21. Extra Tank
- 22. Tissue Basket

Package Info:

Each instrument packed Wooden Case

Outer Case Dimensions: 1300*540*560mm,

Gross Weight: 80 kg