

MTP-1500 Automatic Tissue Processor



Feature :

1. English menu, man-machine interaction, simple operation.
2. Prompt judge and automatic processing in case of abnormal condition, which can prevent the tissue from being damaged.
3. Double-control heat without water
4. Based on different tissues, several running programs available for customers
5. Self-provided power supply can ensure to continue 6 hours' working even if the electricity has been cut off.
6. Environment friendly: The air as purified through internal circulation can be exhausted off the window.
7. With transmission of gear and rack, mechanical failures can be reduced.
8. Several standby programs available

Technical data:

9. Tanks : 12 Tanks (3 wax Tanks in the end)
10. Capacity of single tank : 1500ml
11. Basket can hold 75 specimens
12. Single-tanks processing time: Randomly set within the range from 0 to 99 hour(s).

13. The temperature of wax cylinder can be randomly set within the range of 99□.
14. Time for cylinder to be lifted: 38 seconds
15. Tissue-protection cylinder sequence: 7th cylinder
16. Times of mixing per minute: 3 times

Accessories:

Standard Setting List

17. 1. Tanks 12 pcs
18. 2. Tissue basket 1 pc
19. 3.The exhaust pipe 1 pc
20. 4.Tissue cassette 20 pcs
21. 5. Fuse (5 A) 2 pcs
22. 6. Power line 1 line
23. Optional Accessories
24. Extra Tank
25. Tissue Basket

Package Info:

Each instrument packed Wooden Case

Ourter Case Dimensions: 1310*540*580 mm,

Gross Weight: 100 kg