

## MTP-2000 Automatic Tissue Processor



### Product Features:

- 1) With 16 sets of editable programs, intelligent operation and operation.
- 2) Power failure protection function, automatically switch to the built-in power supply after power failure, up to 8 hours, also configure high-power battery, to ensure that the program is completed in one time.
- 3) Add specimens at any time during the running of the program; you can switch the paraffin bath to warm the wax at any time. Real-time display of program content, heating status, set temperature, real-time temperature and other operating status.
- 4) The scheduled boot time can be set at any time according to the needs of the work.
- 5) The order of the baskets can be set according to actual needs.
- 6) The surface of the machine body adopts APS static electricity, spraying process, high strength, corrosion resistance, streamlined appearance and easy to clean.
- 7) Fully sealed dehydration operation, external exhaust device ensures a green and environmentally friendly working environment.
- 8) There are three sets of temperature control systems to ensure that one set is in normal use and two sets are on standby. Real-time display of temperature during operation Heating status at a glance, temperature control accuracy  $\pm 1^{\circ}\text{C}$
- 9) Automatic switch machine can set the time arbitrarily within one year

### Technical Parameters:

- 1) 12 liquid cylinders (3 of which are wax cylinders)
  - 2) Single cylinder volume: 2200ml
  - 3) Drip retention time: adjustable from 5 to 999 seconds
  - 4) Single cylinder working time: arbitrarily set within 0-24 hours
  - 5) Temperature setting: 0-99 °C arbitrary setting
  - 6) Starting time of the cylinder is 45 seconds - 120 seconds
  - 7) Organization protection cylinder position 5 cylinder-7 cylinder
  - 8) Mixing times: 3 times / minute
  - 9) Power supply: AC220V±10%, 50HZ±10%
  - 10) Power: ≤500W
- Product Function:

**Package Info:**

- © Each instrument packed Wooden Case
- © Outer Case Dimensions: 1250\*410\*560mm,
- © Gross Weight: 128 kg