

STP-500 Spin Tissue Processor



Features:

1. Adopting computer to control whole work process, operate brief work stabilization.
2. Having exchange jar and change place control draw perch, operate agility brief.
3. If lose electric give an alarm, protecting tissue in liquid, obstache, protect function and so on.
4. Adopting agitation dehydrate mode, so tissue and solvent, paraffin meet, thereby well effect.
5. Screen showing information is English, easy to use.
6. Time design range wide: delay time 59hr59min. Each jar working hours longest 59hr59min, Time design least distance 1 min.
7. Wax jar adopt constant temperature control system. Temperature precision high: defend pollution, defend canker.
8. Adopting 9 slices 1.2L medicine aquarium, at any moment observe tissue change at work.
9. According to user's need to tailor standby electrical source, after power cut it can work in great, working hours can reach 4—16hours.

Technical Data:

10. Steps of processing tissue: 12steps.
11. Glass container: Quantity—9capacity-1.2L.
12. Wax cup: quantity—3capacity—1L.

13. Temperature range—45~85 \square \pm 3 \square .
14. Process: Each cup working hours—0~59hr59min.
15. Delay process—0~59hr59min.
16. Vertically reciprocative motion of tissue basket: 3 times/min
17. Tissue basket specification: Barrel tissue basket-- ϕ 95x80mm.
18. Layer nacelle--- ϕ 95x80mm Five layers.
19. Weight: <60kg.
20. External dimension: (diameter)660 x (High)550 x 680(Tiptop)mm.
21. Rated voltage: 220v,50Hz.

Accessories:

Standard Setting List

22. Reagent bottle 9 pcs
23. Paraffin wax tank 3 pcs
24. Paraffin wax tank connection line 3 pcs
25. Tissue basket 1 pc
26. Tissue cassette 20 pcs
27. 6. Hexagonal Wrench 1 pc
28. 7. Screwdriver 1 pc
29. 8. Fuse (5A) 2 pcs
30. 9. Power line 1 line

Optional Accessories

- Extra paraffin tank
- Extra reagent tank

Package Info

Each instrument packed Wooden Case

Outer Case Dimensions (for machine): 670*550*680 mm,
Outer Case Dimensions (for tanks) : 510*510*510mm,

Gross Weight: 110 kg