
MES-700L Tissue Embedding System



Features :

1. Automatic program control, It can be setting the operation time for any day, any time, any minute
2. in a week
3. .The machine is composed of two molds,main console and cold plate. It is independent design,can put the position left to-right or right to -left
4. The large 6L capacity allows user processing 400 tissues at the same time.
5. With manual operation and pedal switch for the operator to choose
6. LED light for embedding space
7. Tweezers storage cylinder is convenient for user
8. Touching LCD scree is simple and convenient, and it is clear at a glance

Technical Data :

9. Main Console:

| | |
|----------------------|-----------------------|
| Model No.: | 10. RD-99 |
| Wax cylinder volume: | 11. 6L |
| Dimensions: | 12. 600 x 650 x 450mm |
| Weight: | 13. 40kgs. |

| | |
|--|--------------------------|
| Voltage: | 14. 110v-240v; 50Hz-60Hz |
| Power: | 15. 1500w |
| Heater temperature setting range: | 16. ambient to 99 ° |
| Temperature for the heat preservation tank | 17. ambient to 99 ° |
| Temperature for the operation table: | 18. ambient to 99 ° |
| Working environmental temperature: | 19. 0 ° - 45 ° |

20. Cryo Console

| | |
|--|--------------------------|
| The size for the cooling surface: | 21. 365x330mm |
| The size for the whole unit: | 22. 380x650x450mm |
| Voltage: | 23. 110v-240v; 50Hz-60Hz |
| Power: | 24. 120W |
| The temperature for the cooling plate: | 25. -20 degree |

Accessories:

Standard Setting List

1. Foot Switch 1 pc
2. Fuse (5 A) 2 pcs
3. Power line 2 line

Optional Accessories

Temperature sensor

Package Info:

Each instrument packed Wooden Case

Ourter Case : 1180*750*620 mm,

Gross Weight : 100 kg