


Immuno Autostainer MIS-3500



1. Open system, Suitable for various kinds of specimen
2. True hands free, walk away ,overnight operation.
3. Easy Windows™ based programming, password protection
4. Allows for several strainers using single computer simultaneously
5. User-specified protocols/User's own reagents, serval runs per day
6. Robust, proven construction and operation
7. No disposable pipette tips
8. Per-run reagent volume checks
9. Detailed reports stored
10. No additional consumable costs
11. No special reagents-using for different company

No	Specification	Model
1	Immuno Autostainer	HS2300
2	Slides Capacity	48
3	Reagents Capacity	64
4	Minimum dropping reagent	50ul
5	Dropping reagent increment	50ul, 100ul, 150ul, 200ul

6	Maximal output power:	≤350VA
7	Working noise	≤60dB(A)
8	Application	IHC,HE more
9	Dimension(mm)	1015*665*520
10	Power supply	220V or 110V
11	Rated current	1.6A
12	Weight	116 KGS