

MSS-600 Slide Stainer

Tissue dyeing process is essential in doing the pathological research, this instrument is designed for saving the fussy process of manual tissue dyeing



Features:

Tissue dyeing process is essential in doing the pathological research, this instrument is designed for saving the fussy process of manual tissue dyeing

Features:

1. High-definition & sensitive of touch LCD
2. Unique design for dyeing process
3. Flexible transmission, noiseless
4. Intellectualized inflow, outflow & drainage system
5. Special dyeing shelf loading with decades of specimen
6. Streaming water cleaning & tissue roasting function

Technica Data:

7. Reagent Containers: 12/16
8. One cleaning bath (the 9th one)
9. One drying bath (the 1st one)
10. Processing time of each container: 0s~59 min and 59 s.
11. Capacity of containers: 1 L, specimen process: 72 pc.
12. Working voltage: AC220 V 50Hz ,AC110 V 60Hz
13. Power: 1000 VA

14. Dimensions: 1130 x 460 x 400 mm

15. Weight: 85 kg

Accessories:

Standard Setting List

16. 1. Tanks 14/18 pcs

17. 2. Slide basket 1 pc

18. 3. Tap 1 set

19. 4. Fuse (5 A) 2 pcs

20. 5. Power line 1 line

21. Optional Accessories

22. Extra Tank

23. Tissue Basket

Package Info:

Each instrument packed Wooden Case

Outer Case Dimensions: 1350*450*400 mm,

Gross Weight: 110 kg