

MSS-800 Slide Stainer



Performance characteristics

1. The surface of the machine body adopts APS static electricity, spraying process, high strength, corrosion resistance, streamlined appearance and easy to clean.
2. It adopts the combination of wire rope transmission and computer control technology. The system runs stably, the positioning is accurate, and the noise is small.
3. The system can store 99 sets of different programs, each with 50 sets of program content, which can meet the special requirements of users for dyeing different tissues.
4. The processing time of a single cylinder can be set arbitrarily between 0 and 99 minutes and 59 seconds.
- 5, dyeing cylinders can be insulated; effectively overcome the shortcomings of poor dyeing caused by low ambient temperature.
- 6, cleaning tank (two); double washing tank is used to strengthen the cleaning of hematoxylin; make the residual hematoxylin cleaning more clean;
7. Automatic cleaning device can be added when changing each cylinder of the reagent cylinder to make the operation more convenient and quick.
8. With air discharge device to reduce air pollution.
- 9, change the tank drain: cleaning reagent time 0-99 seconds

The main technical parameters

- 1, the number of liquid cylinders: 19 (double washing, baking cylinders);
- 2, single cylinder volume: 177mm × 39mm × 98mm;
- 3, the number of single dyeing: 52 pieces / time;

- 4, the number of reagent cylinder insulation: 16;
- 5, reagent tank insulation settings: room temperature -40 ° C;
- 6, single cylinder processing time: can be arbitrarily set within 0 ~ 59 minutes and 59 seconds;
- 7, change the cylinder to stay drain time: can be set arbitrarily within 0 ~ 99 seconds;
- 8, the power supply voltage: AC220V ± 10%;
Frequency: 50±1HZ;Power: less than 1000W;
- 9.Size: 143cm × 44cm × 50cm

Accessories:

Standard Setting List

- Host 1pc
- Dyeing rack 1pc
- Outlet pipe (2M) 1pc
- Inlet pipe (dedicated) 1pc
- Outlet pipe tight circle 1pc
- Exhaust hose 1pc
- Insurance tube (5A) 2pcs
- Environmental protection cover 1pc
- Power cord 1pc
- Manual 1pc
- Instrument and equipment Manual 1pc